

## How to utilize EVA:

On boot, the user will be prompted to enter the wifi network name password and then enter the user's name via the onscreen touch keyboard. From there it will enter the home screen. The user should scan their prescription bottles into the machine and assure that the extracted data is accurate. Any incorrect values can be fixed by touching the incorrect value, entering the correct value using the onscreen keyboard, then hitting enter, then update. Complete these steps for each of your prescriptions.

Once prescriptions are scanned, the device will update the user at the appropriate time to take their medicine. Upon doing so, the user will need to confirm with the device if they have taken their medicine yet. They can do this by selecting: "Yes, No IDK" on the automatically generated GUI. The answer is logged in the data base.

If the user wants drug information, they can tap the drug info button, and then tap the specific drug they want information on. Upon doing so, EVA will scan drugs.com for info and display it to the user.

The Instructions are very simple since this is geared toward elderly individuals. There instructions encapsulate the full functionality of the device for the 1.0 version of EVA.